

# **BSI sCMOS Camera**

# Dhyana95 V2

# Smaller, but Mightier.



#### Smaller size but more functions

Dhyana95 V2 adopts Tucsen's new technology and performance standard. The smallest dimension among peers for scientific imaging makes it more popular for demanding compact space. The water-cooling method and CameraLink interface apply the camera for conditions requiring high stability.

	Dhyana95	Dhyana95 V2
Appearance		O
Air-cooling	•	•
Water-cooling		•
USB3.0	•	•
CameraLink		•
Size (mm)	120x119x121	100x118x127


#### Faster readout by doubling speed

Dhyana95 V2 adds a STD high speed readout mode, of which the frame rate is up to 48fps@4.2MP which is twice as the normal mode. It can be achieved grogressively by using ROI function for applications demanding special frame rate.

ROI	Normal(16-bit)	STD(12-bit)
2048 x 2048	24	48
1608 x 1608	31	61
1200 x 1200	41	81
1024 x 1024	48	95
2048 x 512	95	190
2048 x 256	189	375
2048 x 128	369	739


#### Superior background with accurate calibration


Dhyana95 V2 now can provide a superior background for quantitative analysis applications. The DSNU/PRNU have reached the international optimum of 0.2e- and 0.3% respectively, as the problems from sensor production process, such as edge glow, dead pixels as well as pixel non-uniformity, all have been calibrated accurately.



# 95% quantum efficiency, and excellent SNR


The Dhyana95 V2 uses backside-illuminated sCMOS thinned chip technology to avoid light interference from the wiring layer. It has excellent response capabilities in the ultraviolet, visible and near-infrared bands, with a peak quantum efficiency of up to 95%@560nm, and when the number of incident photons are more than 3 , it can have a better SNR performance than the typical EMCCD camera of 97% peak QE .






### 100Ke- full-well capacity, high dynamic range

Thanks to the 11um large pixel, Dhyana95V2's full-well capacity can be up to 100,000 electrons. Advanced High Dynamic (HDR) mode, through dual-channel gain conbined, it can present high-quality images with rich details of dark and bright in real time, which is very useful in applications such as life science.



# 2" area, large field of view

The 2" array can not only adapt to more optical interfaces and deliver a greater field of view, but also results in fewer lens switches to find the area of interest on the sample.








#### STORM Imaging

"Extensive tests and comparisons with other top-of-the-line EMCCD and sCMOS cameras have been performed in our laboratory. We found Dhyana95, the new sCMOS camera stood its own ground remarkably well and offered satisfactory performance across the board".

—Professor Ning Fang at Georgia State University



#### Soft X-ray Detection

The distinctive features of Dhyana95 make it an excellent alternative to back-illuminated CCDs which are commonly applied in soft X-ray applications of previous-generation synchrotrons, such as coherent scatter experiments. They would greatly benefit fully from the Dhyana95's high frame rates.

——Kewin Desjardins, from the French SOLEIL synchrotron





Figure 1. A 5 micron pinhole diffraction spot of 186 eV recorded by Dhyana95, which shows excellent dynamic range with its image diffraction order of the reached to the maximum value, the 6th.

Figure 2. A composited radia-illumination image accumulated from 50 images and each with a 100 ms exposure time by the Dhyana95. The total acquisition time of 50 images is less than 10 s, while the same image effect achieved with a backilluminated CCD camera is a couple of minutes.

Model	Dhyana95 V2
Sensor size	2"
Sensor model	Backside-illuminated sCMOS
Color/Monochrome	Monochrome
Quantum efficency	95%@560nm
Effective pixels	2048(H) x 2048(V)
Pixel size	11um x 11um
Effective area	22.5mmx22.5mm
Full well capacity	80000e-@HDR,100000e-@STD(Low gain)
Frame rate	24fps@16bit HDR,48fps@12bit STD
Read noise	1.6e-(median), 1.7e-(rms)
Shutter type	Rolling Shutter
Exposure time	21us-10s
DSNU	0.2e-
PRNU	0.3%
Offset	100 (HDR and STD)
Cooling method	Forced air: (Ambient at +20°C):-15°C Water: (Ambient at +20°C): -25°C
Dark current	Forced air: 0.5e-/pixel/s, Water: 0. 25e-/pixel/s
Binning	1x1, 2x2, 4x4
Sub-array	2048x1024, 2048x512, 1608x1608, 1200x1200 1024x1024, 512x512, 256x256
External trigger mode	Standard/Synchronous/Globaltrigger/Software trigger
Trigger delay function	0-10s(1µs steps)
Trigger output	Exposure/Global/readout signal/High/Low
External trigger routing	SMA
Timestamp Accuracy	1us
Digital interface	
J	USB3.0, CameraLink
SDK	USB3.0, CameraLink C / C++
SDK	C/C++
SDK Bit depth	C / C++ 16bit &12bit
SDK Bit depth Lens mount	C / C++  16bit &12bit  C-mount & F-mount
SDK Bit depth Lens mount Power supply	C / C++  16bit &12bit  C-mount & F-mount  12V/8A
SDK  Bit depth  Lens mount  Power supply  Power consumption	C / C++  16bit &12bit  C-mount & F-mount  12V/8A  ~60W  C-mount 100mmx118mmx127mm
SDK Bit depth Lens mount Power supply Power consumption Camera size	C / C++  16bit &12bit  C-mount & F-mount  12V/8A  ~60W  C-mount 100mmx118mmx127mm F-mount 100mmx118mmx157mm  Mosaic 1.6 / labview / Matlab / metamorph /
SDK  Bit depth  Lens mount  Power supply  Power consumption  Camera size  Software	C / C++  16bit &12bit  C-mount & F-mount  12V/8A  ~60W  C-mount 100mmx118mmx127mm F-mount 100mmx118mmx157mm  Mosaic 1.6 / labview / Matlab / metamorph / Micromanager

#### **Dimensions**



#### **Packing list**

Camera * 1	USB3.0 cable * 1
Power cord * 1	F optical interface * 1
Power adapter * 1	Software (U disk) * 1
Water-cooled pipe * 2 (opt.)	Hexagon wrench *1
CameraLink cable * 2 (opt.)	

#### **Tucsen Photonics Co., Ltd.**

Address: 5#, No.2 Yangqi Branch Rd, Gaishan Town, Cangshan Area, Fuzhou, Fujian,PRC,China.

Tel: + 86-591-28055080 Website: www.tucsen.com Email: support@tucsen.com

